

Our priority is your total satisfaction. Should you ever have a problem with any standard, Water QC or otherwise, let us know. We'll immediately investigate the problem by testing a retained sample of your solution. If the error is on our end, you'll be offered a full refund or a free replacement — your choice.



Contents

Water Standards	100
Need a Custom CRM?	16

- ✓ **Five-year shelf life**
- ✓ **Traceable to NIST SRMs**
- ✓ **Produced under ISO 9001**
- ✓ **Produced under ISO 17025**
- ✓ **Produced under ISO 17034**
- ✓ **Assayed by optimal validated procedures**

Custom potable water standards for certain products are available upon request.

Bromate	
ICBR031	Matrix: H ₂ O
ICBR031-125ML ICBR031-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
BrO ₃ ⁻	1,000

Chlorate	
ICCL031	Matrix: H ₂ O
ICCL031-125ML ICCL031-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
ClO ₃ ⁻	1,000

Bromide	
ICBR1	Matrix: H ₂ O
ICBR1-125ML ICBR1-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
Br ⁻	1,000

Chlorite	
ICCL021	Matrix: H ₂ O
ICCL021-125ML ICCL021-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
ClO ₂ ⁻	1,000

NOTE: Contains less than 10ppm ClO₃⁻.

1,000 µg/mL Total Cyanide	
CN-1000-25-20ML	Volume: 20 mL Matrix: H ₂ O
Analyte	µg/mL
CN ⁻	1,000

For the determination of cyanide in aqueous samples.

Mercury Standard	
MSHG-1PPM	Matrix: HCl
MSHG-1PPM-125ML MSHG-1PPM-500ML	Volume: 125 mL Volume: 500 mL
Analyte	µg/mL
Hg	1

Custom wastewater standards are available upon request.

ANALYTE	MATRIX	STARTING MATERIAL	µg/ML	VOLUME	CATALOG #
Total Organic Carbon, TOC	H ₂ O	KHP	1,000	125 mL 500 mL	TOCKHP1-125ML TOCKHP1-500ML